EN

## Mediblink Spacer for Aerosol Medication MediFlow

### INSTRUCTIONS FOR USE

#### 1. Important Information

##### 1.1 General

Please read these instructions for use carefully and completely. Do not discard them, so you can consult them at a later date. If you fail to comply with the instructions for use, injury or damage to the device cannot be ruled out.

**If symptoms persist or if your condition worsens, discontinue the application and contact your doctor immediately.**

##### 1.2 Structure of safety instructions

Safety-critical warnings are categorized according to hazard levels in these instructions for use:

- The signal word WARNING is used to indicate hazards which, without precautionary measures, can result in serious injury or even death.
- The signal word CAUTION is used to indicate hazards which, without precautionary measures, can result in minor to moderate injury or impair treatment.
- The signal word NOTICE is used to indicate general precautionary measures which are to be observed to avoid damaging the product during use.

##### 1.3 Treatment of anyone who requires assistance

Anyone who requires assistance must be supervised constantly by an adult during inhalation therapy.

This is the only way to ensure safe and effective treatment.

Individuals in this group often underestimate the hazards involved, thus resulting in a risk of injury.

The product contains small parts. Small parts can block the air-ways and lead to a choking hazard. Therefore, make sure that you always keep all components of the product out of the reach of babies, infants and children.

This product is only designed for patients who are able to breathe by themselves and are conscious. This is the only way to ensure effective treatment and avoid a choking hazard.

##### 1.4 Hygiene

Only use the MediFlow spacer and accessories that have been cleaned and dried as described in chapter 4 for inhalation therapy.

Contamination and residual moisture encourage the growth of bacteria, so increasing the risk of infection.

Therefore, please observe the following hygiene instructions:

- Each MediFlow spacer must only be used at home by a single patient for hygiene reasons.
- Make absolutely sure you also carry out cleaning and drying before using the device for the first time.
- For cleaning and disinfection, always use fresh tap water of drinking water quality.
- Make sure all components are dried properly whenever you have performed cleaning or disinfection.
- Do not keep the MediFlow spacer and accessories in a damp environment or together with damp objects.

#### 2. Product Description

##### 2.1 Components

Check that all components of your MediFlow spacer product are contained in your package. If anything is missing, please notify the dealer from whom you purchased the product immediately.

| No. | Components             | Specification/Size    |                        |                       |                |
|-----|------------------------|-----------------------|------------------------|-----------------------|----------------|
|     |                        | L<br>(M493–<br>Large) | M<br>(M492–<br>Medium) | S<br>(M491–<br>Small) | N<br>(No mask) |
|     |                        | ✓                     | ×                      | ×                     | ×              |
| 1.  | Masks                  | ×                     | ✓                      | ×                     | ×              |
|     |                        | ×                     | ×                      | ✓                     | ×              |
| 2.  | Dust cap<br>(optional) | ✓                     | ✓                      | ✓                     | ✓              |
| 3.  | Mouthpiece             | ✓                     | ✓                      | ✓                     | ✓              |
| 4.  | Storage body           | ✓                     | ✓                      | ✓                     | ✓              |
| 5.  | Back piece             | ✓                     | ✓                      | ✓                     | ✓              |
| 6.  | Valves                 | ✓                     | ✓                      | ✓                     | ✓              |
| 7.  | Whistle<br>(optional)  | ✓                     | ✓                      | ✓                     | ✓              |

Note: "✓" stands for having these components while "×" for none.

##### 2.2 Intended purpose

The Spacer for aerosol Medication MediFlow is a holding chamber, and is used in conjunction with medication sprays or "metered dose inhalers" (MDIs) in the treatment of diseases of the airways.

##### 2.3 Intended use

MediFlow spacer is intended for use by treating chronic obstructive pulmonary disease (COPD) and asthma with medicine-assisted inhalation device.

##### 2.4 Combination of medical devices

The MediFlow spacer should be used in conjunction with "metered dose inhalers" (MDIs) in the treatment of diseases of the airways. The metered dose inhalers are approximate rectangular, 23±2mm long x16.5±2mm wide.

##### 2.5 Contraindications

There are no contraindications known.

##### 2.6 Product combinations

The MediFlow spacer has four specifications: 'S' is Small mask (Reference number M491/ model DL-01D-I), 'M' is Medium mask (Reference number M493/ model DL-01D-C), 'L' is Large mask (Reference number M492/ model DL-01D-L), and N is without mask. Patients can choose the appropriate specifications based on the size of the face.

##### 2.7 Description of function

The MediFlow spacer helps to minimize coordination mistakes when using metered dose inhalers, avoid deposition of large quantities of medication in the oropharyngeal area, and the undesirable side effects associated there with a flexible connecting adaptor. It can be used with most conventional metered dose inhalers.

##### 2.8 Material information

| Materials            | Parts                          |
|----------------------|--------------------------------|
| PETG/PP/ABS          | Inhalation chamber, mouthpiece |
| SILICONE             | Valve in the chamber, mask     |
| THERMOPLASTIC RUBBER | Back piece, dust cap           |

All materials are latex-free and BPA-free.

##### 2.9 Shelf life

The shelf life of the MediFlow Spacer for Aerosol is 3 years. The product is recommended to be replaced after 12 months or to be replaced after using it for 1460 times.

#### 3. Inhalation

##### 3.1 Function check

###### Clear pathway through the valve

**After receiving the MediFlow spacer**, check that the path through the valve in the mouthpiece is unobstructed.

If the valve is obstructed, the MediFlow spacer must not be used.

###### Position of the valve

**Before every use**, check that the valve is in the correct position inside the MediFlow spacer mouthpiece:

- Check that the valve is in the position shown inside the mouthpiece. If necessary, use a stream of water to correct its position.
- Dry the MediFlow spacer completely before using it again.

##### 3.2 Preparation for inhalation

###### WARNING

Since the inhalation chamber of the MediFlow spacer is not completely closed, small particles can get in, and might be breathed in during an inhalation session (choking hazard). Therefore, it is imperative to make sure there are no foreign bodies inside the MediFlow spacer before every use.

###### CAUTION


Check all product components and the accessories before each use. Replace any broken, misshapen or seriously discolored parts. Please also follow the instructions for assembly given below. Damaged components and an in-correctly assembled holding chamber may impair the function of the holding chamber and thus the treatment as well.

###### CAUTION


If a mask is being used, make sure that the expiratory valve shim is pressed out. If the valve shim is deflected inward, the dosage delivered may be insufficient when inhaling.

##### 3.3 Performing the inhalation

1) Carefully examine the MediFlow spacer for damage, missing parts or foreign objects. Remove any foreign objects prior to use.

2) Remove cap from the chamber (if applicable) and place the mask connector over chamber's mouthpiece until it fits securely (if applicable).


###### CAUTION

Before the patient starts inhaling, ensure that all parts are firmly connected to each other. Otherwise, insufficient dosage may be delivered.

3) Remove cap from the Meter-dose inhaler (MDI), shake MDI immediately before use and insert MDI into the backpiece of the chamber.

4) Try to exhale, and hold mask against the face, press MDI and inhale slowly. Hold mask in place and breath in and out through the chamber for 5-6 breaths. If there is no mask, place lips around mouthpiece, press MDI and inhale slowly. Slow down inhalation if a tone is heard. Hold breath for 5-10 seconds.

###### CAUTION

If you use a mask, make sure that the mask completely covers both corners of the mouth and the nose. Otherwise, the treatment can be less effective due to escaping aerosol, because of underdosage, for example. For possible side effects of escaping aerosol, please read the instructions for use of the medication in question.

###### CAUTION

Children younger than 18 months should inhale four to six times; children older than 18 months should inhale deeply two to four times. Otherwise, insufficient dosage may be delivered. For adults, one breath is enough to inhale the medication from the MediFlow spacer.

5) If the doctor asks you to take medicine for more than one spray, please take a rest for one minute, and then repeat steps 3-4.

#### 4. Cleaning


Regularly clean all product components and the accessories used, if they are visibly dirty immediately after use. The chamber is usually cleaned once a week.

##### 4.1 Preparation

- Take the metered dose inhaler (MDI) out of the connecting ring on the inhalation chamber and close it with the protective cap provided for this purpose.
- If applicable, detach the mask from the mouthpiece.
- Disconnect the mouthpiece from the inhalation chamber (if applicable).

##### 4.2 Cleaning


- Remove the back piece and disassemble the other parts of the MediFlow spacer.

- Soak all parts for 15 minutes in lukewarm water (temperature must not exceed 60°C) with liquid detergent. Agitate gently.

- Rinse in clean water.


- Shake out excess water. Do not rub dry.

- Let air dry in vertical position.

- Re-assemble back-piece when unit is completely dry and ready for use.

NOTE: If you wash the product in the dishwasher, make sure that the temperature does not exceed 60°C.

###### CAUTION

**Do not clean with boiling or hot water as this will damage your MediFlow spacer. Only silicone mask may be cleaned in boiling or hot water.**

##### 4.3 Visual inspection

Inspect all product components after each cleaning and disinfection. Replace any broken, misshapen or seriously dis-colored parts.

##### 4.4 Drying and storage

- Dry all individual parts completely. Because of the antistatic properties of the inhalation chamber, don't rub dry or heat dry.
- Place the protective cap over the mouthpiece (if applicable). Info: The protective cap protects the mouthpiece from damage during transport and storage.
- Attach the mouthpiece to the inhalation chamber (if applicable).
- Store the MediFlow spacer and its accessories in a dry, dust-free place (safely out of continuous direct sunlight).

#### 5. Miscellaneous

##### 5.1 Disposal

All product components can be disposed of with domestic waste unless this is prohibited by the disposal regulations prevailing in the respective member countries.

##### 5.2 Explanation of symbols

The following symbols are displayed on the packaging of your product MediFlow spacer:

| Symbol | Explanation                     |
|--------|---------------------------------|
|        | Manufacturer                    |
|        | Authorized representative in EU |
|        | Importer                        |

| Symbol | Explanation                                                                                                                                                                                                                                                                                                                           |
|--------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|        | Distributor                                                                                                                                                                                                                                                                                                                           |
|        | CE Marking of Conformity, Notified body ID number                                                                                                                                                                                                                                                                                     |
|        | Before use, read the instructions. Electronic instructions for use: <a href="http://www.mediblink.com/t/m491.pdf">http://www.mediblink.com/t/m491.pdf</a> <a href="http://www.mediblink.com/t/m492.pdf">http://www.mediblink.com/t/m492.pdf</a> <a href="http://www.mediblink.com/t/m493.pdf">http://www.mediblink.com/t/m493.pdf</a> |
|        | Warning                                                                                                                                                                                                                                                                                                                               |
|        | Reference number                                                                                                                                                                                                                                                                                                                      |
|        | Batch number                                                                                                                                                                                                                                                                                                                          |
|        | Use-by date                                                                                                                                                                                                                                                                                                                           |
|        | Model number                                                                                                                                                                                                                                                                                                                          |
|        | Medical device                                                                                                                                                                                                                                                                                                                        |
|        | Unique Device Identifier                                                                                                                                                                                                                                                                                                              |
|        | Number of products included in packaging                                                                                                                                                                                                                                                                                              |
|        | Paper                                                                                                                                                                                                                                                                                                                                 |
|        | Keep dry                                                                                                                                                                                                                                                                                                                              |

\*The production date can be read from LOT number [YYYYMM]; first four digits represent the Year and last two digits Month of production. Example: LOT 202503 = March 2025

##### 5.3 Contact

For all product information and in the event of defects or questions about usage, please contact your local Mediblink distributor.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

**Manufactured for (Importer for EU & Distributor):**  
Mediblink d.o.o., Gubčeva cesta 19, 8210 Trebnje, Slovenia, info@mediblink.com  
www.mediblink.com

**Tai'an Dalu Medical Instrument Co., Ltd.**  
No.3988 Yitianmen Street, Hi-tech Zone, Tai'an, 271000 Shandong, P.R. China  
E-mail: sales@taiandalumed.com  
**Lotus NL B.V.**, Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands.  
E-mail: peter@lotusnl.com

Version: A/6  
Issue Date: 2.12.2024  
Date of last revision: 25.2.2025  
1.0, 25.2.2025

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